



Vital Signs AGC

API Documentation
revision v2.2.0

Generated by Doxygen

Contents

1	Vital Signs AGC	1
1.1	Overview	1
2	Module Documentation	3
2.1	AGC Module	3
2.1.1	Detailed Description	4
2.1.2	Macro Definition Documentation	5
2.1.3	Typedef Documentation	5
2.1.4	Enumeration Type Documentation	5
2.1.5	Function Documentation	6
2.2	HAL Functions	11
2.2.1	Detailed Description	11
2.2.2	Function Documentation	11
2.3	Error Codes	16
2.3.1	Detailed Description	17
2.3.2	Typedef Documentation	17
2.3.3	Enumeration Type Documentation	17
3	Data Structure Documentation	19
3.1	agc_configuration_t Struct Reference	19
3.1.1	Detailed Description	19
3.1.2	Field Documentation	19
3.2	agc_status_t Struct Reference	21
3.2.1	Detailed Description	21
3.2.2	Field Documentation	21
3.3	agc_version_t Struct Reference	21
3.3.1	Detailed Description	22
3.3.2	Field Documentation	22

1 Vital Signs AGC

1.1 Overview

The Automatic Gain Compensation (AGC) algorithm keeps the PPG signal within a configured ADC range with optional amplitude control. The PPG signal range is controlled via the PD offset and the PPG amplitude is controlled via the LED current.





Figure 1 Block Diagram

PD offset control is self-adjusting regarding the impact of PD offset changes and ADC sensitivity. The integrated LED current control algorithm follows a simple statistical approach to identify phases of a steady PPG signal, free of heavy motion artifacts. A proper configuration is required for good algorithm performance. The module also supports that the LED current is controlled externally by a Vital Signs algorithm.

The AGC algorithm is a single PPG channel solution that can be applied independently to multiple PPG channels, if the sensor supports setting PD offset and LED current per PPG channel. This means that this solution is suitable for single-channel HRM, multi-channel HRM, and dual-channel SpO2 applications when used with a compatible sensor such as AS7050 and AS7056/57. AS703x sensors are not suited for multi-channel applications with this AGC algorithm.

**Figure 2 State Diagram**

On a high-level overview, the AGC module is used as follows:

1. After initialization, the AGC module enters Configuration state and can be configured based on the sensor configuration and the application.
2. After the configuration has been set, the module enters the Processing state. When entering the Processing state, the HAL layer and the internal parameters of the module are initialized and the AGC algorithm is ready for execution.
3. The function to execute the algorithm can be called anytime while in the Processing state. The algorithm processes the provided input data and generates the resulting PD offset and LED current values. The resulting PD offset and LED current values are applied to the chip via HAL functions.
4. The AGC status can be read anytime while in the Processing state.
5. When calling the function to stop processing, the AGC module enters the Configuration state. Internally, this function de-initializes the HAL layer.
6. The configuration of the AGC module can be updated while in the Configuration state.




Figure 3 Sequence Diagram

2 Module Documentation

2.1 AGC Module

Data Structures

- struct `agc_status_t`
- struct `agc_configuration_t`
- struct `agc_version_t`

Macros

- `#define AGC_MAX_CHANNEL_CNT 4`

Typedefs

- `typedef uint8_t agc_change_state_t`
- `typedef uint8_t agc_mode_t`
- `typedef uint8_t agc_ampl_cntl_mode_t`
- `typedef uint8_t agc_channel_id_t`

Enumerations

- `enum agc_change_state {`
 `AGC_STATE_UNCHANGED = 0,`
 `AGC_STATE_INCREASED = 1,`
 `AGC_STATE_DECREASED = 2,`
 `AGC_STATE_ATMIN = 3,`
 `AGC_STATE_ATMAX = 4,`
 `AGC_STATE_NOT_CONTROLLED = 5,`
 `AGC_STATE_UNDEFINED = 6 }`
- `enum agc_mode {`
 `AGC_MODE_DEFAULT = 0,`
 `AGC_MODE_STATUS = 1 }`
- `enum agc_ampl_cntl_mode {`
 `AGC_AMPL_CNTL_MODE_DISABLED = 0,`
 `AGC_AMPL_CNTL_MODE_AUTO = 1,`
 `AGC_AMPL_CNTL_MODE_EXTERNAL = 2 }`

Functions

- `err_code_t agc_initialize (void)`
- `err_code_t agc_set_configuration (const agc_configuration_t *p_agc_configs, uint8_t agc_config_num)`
- `err_code_t agc_get_configuration (agc_configuration_t *p_agc_configs, uint8_t *p_agc_config_num)`
- `err_code_t agc_start_processing (void *p_config, uint16_t size)`
- `err_code_t agc_execute (const uint8_t *p_fifo_data, uint16_t fifo_data_size)`
- `err_code_t agc_get_status (agc_status_t *p_agc_status, uint8_t *p_agc_status_num)`
- `err_code_t agc_stop_processing (void)`
- `err_code_t agc_shutdown (void)`
- `err_code_t agc_get_version (agc_version_t *p_agc_version)`

2.1.1 Detailed Description

The Automatic Gain Compensation (AGC) module implements an algorithm to keep PPG signals within a configured range by continuously adjusting the photodiode (PD) offset current and the LED current.

2.1.2 Macro Definition Documentation

2.1.2.1 AGC_MAX_CHANNEL_CNT `#define AGC_MAX_CHANNEL_CNT 4`

The maximum number of enabled AGC instances.

2.1.3 Typedef Documentation

2.1.3.1 agc_change_state_t `typedef uint8_t agc_change_state_t`

Type for [agc_change_state](#)

2.1.3.2 agc_mode_t `typedef uint8_t agc_mode_t`

Type for [agc_mode](#).

2.1.3.3 agc_ampl_cntl_mode_t `typedef uint8_t agc_ampl_cntl_mode_t`

Type for [agc_ampl_cntl_mode](#).

2.1.3.4 agc_channel_id_t `typedef uint8_t agc_channel_id_t`

Represents a PPG channel. A value of 0 represents a disabled channels, all other values are chip-specific.

2.1.4 Enumeration Type Documentation

2.1.4.1 agc_change_state `enum agc_change_state`

AGC current change information.

Enumerator

AGC_STATE_UNCHANGED	Current not changed.
AGC_STATE_INCREASED	Current increased.
AGC_STATE_DECREASED	Current decreased.
AGC_STATE_ATMIN	Current needs to be decreased, but is at the lower limit.
AGC_STATE_ATMAX	Current needs to be increased, but is at the upper limit.
AGC_STATE_NOT_CONTROLLED	Current is not controlled as amplitude control is disabled.
AGC_STATE_UNDEFINED	Current is undefined. For example the channel has no connected LED.

2.1.4.2 agc_mode enum agc_mode

Operating modes of the AGC algorithm.

Enumerator

AGC_MODE_DEFAULT	Default mode.
AGC_MODE_STATUS	AGC algorithm is disabled but returns at least one status package for the selected channel.

2.1.4.3 agc_ampl_cntl_mode enum agc_ampl_cntl_mode

Amplitude control modes of the AGC algorithm.

Enumerator

AGC_AMPL_CNTL_MODE_DISABLED	No amplitude control enabled.
AGC_AMPL_CNTL_MODE_AUTO	Internal amplitude control enabled.
AGC_AMPL_CNTL_MODE_EXTERNAL	Amplitude is controlled by an external algorithm.

2.1.5 Function Documentation

2.1.5.1 agc_initialize() err_code_t agc_initialize (void)

Initializes the AGC module.

AGC module transitions to Configuration state after initialization.

Return values

ERR_SUCCESS	Function returns without error.
--------------------	---------------------------------

2.1.5.2 agc_set_configuration() err_code_t agc_set_configuration (const agc_configuration_t * p_agc_configs, uint8_t agc_config_num)

Sets the configuration of the AGC module.

All previous AGC configurations are discarded when this function is called. This function can only be called when the AGC module is in Configuration state.

Parameters

in	<i>p_agc_configs</i>	Pointer to an array of AGC configurations. A configuration needs to be provided for each AGC instance that shall be enabled.
in	<i>agc_config_num</i>	Number of AGC configurations in the <i>p_agc_configs</i> array.

Return values

<i>ERR_SUCCESS</i>	Function returns without error.
<i>ERR_ARGUMENT</i>	At least one parameter is wrong.
<i>ERR_POINTER</i>	Detected NULL pointer.
<i>ERR_PERMISSION</i>	Invalid state.

2.1.5.3 agc_get_configuration() `err_code_t agc_get_configuration ( agc_configuration_t * p_agc_configs, uint8_t * p_agc_config_num )`

Gets the configuration of the AGC module.

This function can only be called when the AGC module is not in Uninitialized state.

Parameters

out	<i>p_agc_configs</i>	Pointer to where the AGC configurations are written to. The configuration for each enabled AGC instance is provided.
in, out	<i>p_agc_config_num</i>	Pointer to a variable containing the maximum number of AGC configurations that can be written to the <i>p_agc_configs</i> array. After writing the AGC configurations, the variable is updated to contain the number of written AGC configurations.

Return values

<i>ERR_SUCCESS</i>	Function returns without error.
<i>ERR_SIZE</i>	Destination buffer is too small.
<i>ERR_POINTER</i>	Detected NULL pointer.
<i>ERR_PERMISSION</i>	Invalid state.

2.1.5.4 agc_start_processing() `err_code_t agc_start_processing (`
`void * p_config,`
`uint16_t size )`

Prepares the AGC algorithm for execution.

This function can only be called when the AGC module is in Configuration state. The AGC module transitions to Processing state when this function executes successfully.

Parameters

in	<i>p_config</i>	Chip-specific configuration information, which is forwarded to the HAL.
in	<i>size</i>	Size of the p_config configuration information.

Return values

<i>ERR_SUCCESS</i>	Function returns without error.
<i>ERR_POINTER</i>	Detected NULL pointer.
<i>ERR_SYSTEM_CONFIG</i>	Invalid sample period provided by HAL.
<i>ERR_CONFIG</i>	AGC mode is not correctly configured.
<i>ERR_PERMISSION</i>	Invalid state.

2.1.5.5 agc_execute() `err_code_t agc_execute (`
`const uint8_t * p_fifo_data,`
`uint16_t fifo_data_size )`

Executes the AGC algorithm by processing the provided FIFO data.

This function can only be called when the AGC module is in Processing state.

Parameters

in	<i>p_fifo_data</i>	FIFO data.
in	<i>fifo_data_size</i>	Size of FIFO data.

Return values

<i>ERR_SUCCESS</i>	Function returns without error.
<i>ERR_SIZE</i>	Size of FIFO data is zero.
<i>ERR_POINTER</i>	Detected NULL pointer.
<i>ERR_PERMISSION</i>	Invalid state.

2.1.5.6 agc_get_status() `err_code_t agc_get_status (`
`agc_status_t * p_agc_status,`
`uint8_t * p_agc_status_num )`

Gets the AGC status information for all enabled AGC instances.

This function can only be called when the AGC module is in Processing state.

Parameters

out	<i>p_agc_status</i>	Pointer to where the AGC status informations are written to. The AGC status information for each enabled AGC instance is provided. The order of the AGC status informations is identical to the order used for the configuration of the AGC instances.
in, out	<i>p_agc_status_num</i>	Pointer to a variable containing the maximum number of AGC status informations that can be written to the p_agc_status array. After writing the AGC status informations, the variable is updated to contain the number of written AGC status informations.

Return values

<i>ERR_SUCCESS</i>	Function returns without error.
<i>ERR_SIZE</i>	Destination buffer is too small.
<i>ERR_POINTER</i>	Detected NULL pointer.
<i>ERR_PERMISSION</i>	Invalid state.

2.1.5.7 agc_stop_processing() `err_code_t agc_stop_processing (`
`void )`

Exits the Processing state of the AGC algorithm.

This function can only be called when the AGC module is not in Uninitialized state. The AGC module transitions to Configuration state when this function executes successfully.

Return values

<i>ERR_SUCCESS</i>	Function returns without error.
<i>ERR_PERMISSION</i>	Invalid state.

2.1.5.8 agc_shutdown() `err_code_t agc_shutdown (`
`void )`

De-initializes the AGC module.

Return values

<code>ERR_SUCCESS</code>	Function returns without error.
--	---------------------------------

2.1.5.9 agc_get_version() `err_code_t` agc_get_version (
 `agc_version_t` * *p_agc_version*)

Gets the version of AGC module.

Parameters

out	<i>p_agc_version</i>	Pointer to where the version is written to.
-----	----------------------	---

Return values

<code>ERR_SUCCESS</code>	Function returns without error.
<code>ERR_POINTER</code>	Detected NULL pointer.

2.2 HAL Functions

Functions

- [err_code_t agc_hal_initialize](#) (void *p_config, uint16_t size)
- [err_code_t agc_hal_set_led_current](#) (agc_channel_id_t channel_id, uint8_t led_current)
- [err_code_t agc_hal_get_led_current](#) (agc_channel_id_t channel_id, uint8_t *p_led_current)
- [err_code_t agc_hal_set_pd_offset](#) (agc_channel_id_t channel_id, uint8_t pd_offset)
- [err_code_t agc_hal_get_pd_offset](#) (agc_channel_id_t channel_id, uint8_t *p_pd_offset)
- [err_code_t agc_hal_extract_samples](#) (agc_channel_id_t channel_id, const uint8_t *p_fifo_data, uint16_t fifo_data_size, uint8_t fifo_data_changed, int32_t *p_chan_data, uint16_t *p_chan_data_num)
- [err_code_t agc_hal_get_sampling_period_us](#) (agc_channel_id_t channel_id, uint32_t *p_sampling_period)
- [err_code_t agc_hal_shutdown](#) (void)

2.2.1 Detailed Description

This file consists of platform-dependent interface functions and definitions for AGC module. The HAL needs to be implemented for every chip using the AGC algorithm.

2.2.2 Function Documentation

2.2.2.1 agc_hal_initialize() [err_code_t](#) agc_hal_initialize (
 void * p_config,
 uint16_t size)

Initializes the HAL interface for the AGC.

Parameters

in	<i>p_config</i>	Application-specific configuration.
in	<i>size</i>	Size of the configuration.

Return values

ERR_SUCCESS	Function returns without error.
ERR_SIZE	The size of the configuration does not match.
ERR_POINTER	Detected NULL pointer for data.

2.2.2.2 agc_hal_set_led_current() [err_code_t](#) agc_hal_set_led_current (
 [agc_channel_id_t](#) channel_id,
 uint8_t led_current)

Sets the LED current.

Parameters

in	<i>channel_id</i>	Channel number.
in	<i>led_current</i>	New LED current value.

Return values

<i>ERR_SUCCESS</i>	Function returns without error.
<i>ERR_DATA_TRANSFER</i>	Error during communication with sensor.
<i>ERR_PERMISSION</i>	HAL/OSAL library is not initialized.
<i>ERR_ARGUMENT</i>	LED current is wrong or channel ID is not supported.
<i>ERR_LED_ACCESS</i>	Error on LED configuration or access.

2.2.2.3 agc_hal_get_led_current() `err_code_t agc_hal_get_led_current ( agc_channel_id_t channel_id, uint8_t * p_led_current )`

Gets the LED current.

Note

This function will be called often and it is recommended to cache the LED currents internally to avoid unneeded I2C traffic.

Parameters

in	<i>channel_id</i>	Channel number.
out	<i>p_led_current</i>	Pointer to memory where the LED current will be saved.

Return values

<i>ERR_SUCCESS</i>	Function returns without error.
<i>ERR_DATA_TRANSFER</i>	Error during communication with sensor.
<i>ERR_PERMISSION</i>	HAL/OSAL library is not initialized.
<i>ERR_ARGUMENT</i>	Channel ID is not supported.
<i>ERR_LED_ACCESS</i>	Error on LED configuration or access.
<i>ERR_POINTER</i>	NULL pointer detected for p_led_current.

2.2.2.4 agc_hal_set_pd_offset() `err_code_t agc_hal_set_pd_offset (`
`agc_channel_id_t channel_id,`
`uint8_t pd_offset )`

Sets the PD offset currents.

Parameters

in	<i>channel_id</i>	Channel number.
in	<i>pd_offset</i>	New PD offset value.

Return values

<i>ERR_SUCCESS</i>	Function returns without error.
<i>ERR_DATA_TRANSFER</i>	Error during communication with sensor.
<i>ERR_ARGUMENT</i>	PD offset configuration is wrong or channel ID is not supported.
<i>ERR_PERMISSION</i>	HAL/OSAL library is not initialized.

2.2.2.5 agc_hal_get_pd_offset() `err_code_t agc_hal_get_pd_offset (`
`agc_channel_id_t channel_id,`
`uint8_t * p_pd_offset )`

Gets the PD offset currents.

Parameters

in	<i>channel_id</i>	Channel number.
out	<i>p_pd_offset</i>	Pointer to memory where the PD offset value will be saved.

Return values

<i>ERR_SUCCESS</i>	Function returns without error.
<i>ERR_DATA_TRANSFER</i>	Error during communication with sensor.
<i>ERR_ARGUMENT</i>	Channel ID is not supported.
<i>ERR_PERMISSION</i>	HAL/OSAL library is not initialized.
<i>ERR_POINTER</i>	NULL pointer detected for p_pd_offset.

2.2.2.6 agc_hal_extract_samples() `err_code_t agc_hal_extract_samples (`
`agc_channel_id_t channel_id,`
`const uint8_t * p_fifo_data,`

```
uint16_t fifo_data_size,
uint8_t fifo_data_changed,
int32_t * p_chan_data,
uint16_t * p_chan_data_num )
```

Extract the requested samples from FIFO data stream.

When the AGC module receives a new FIFO data stream, it calls this function for each enabled AGC instance to extract the corresponding samples. When the function is called the first time after new FIFO data has been received, the parameter `fifo_data_changed` is set to TRUE to indicate the FIFO data change.

Parameters

in	<i>channel_id</i>	Channel number.
in	<i>p_fifo_data</i>	FIFO data.
in	<i>fifo_data_size</i>	Size of FIFO data.
in	<i>fifo_data_changed</i>	Set to TRUE to signal that this is new FIFO data.
out	<i>p_chan_data</i>	Pointer to buffer where the sub-sample data can be saved.
in, out	<i>p_chan_data_num</i>	Input: Number of provided <code>p_chan_data</code> elements; Output: Number of used <code>p_chan_data</code> elements.

Return values

<i>ERR_SUCCESS</i>	Function returns without error.
<i>ERR_ARGUMENT</i>	At least one parameter is wrong.
<i>ERR_POINTER</i>	Detected NULL pointer for data.
<i>ERR_PERMISSION</i>	HAL/OSAL library is not initialized.

2.2.2.7 agc_hal_get_sampling_period_us() `err_code_t agc_hal_get_sampling_period_us (`
`agc_channel_id_t channel_id,`
`uint32_t * p_sampling_period )`

Requests the sampling period for the given channel in microseconds.

Parameters

in	<i>channel_id</i>	Channel number.
out	<i>p_sampling_period</i>	Sampling period of the PPG signal in microseconds.

Return values

<i>ERR_SUCCESS</i>	Function returns without error.
<i>ERR_POINTER</i>	If the argument is NULL.
<i>ERR_PERMISSION</i>	HAL/OSAL library is not initialized.

2.2.2.8 agc_hal_shutdown() `err_code_t` agc_hal_shutdown (
void)

Disables this library and block the functions calls.

Return values

<code>ERR_SUCCESS</code>	Function returns without error.
--------------------------	---------------------------------

2.3 Error Codes

Typedefs

- typedef enum `error_codes` `err_code_t`

Enumerations

- enum `error_codes` {
 `ERR_SUCCESS` = 0,
 `ERR_PERMISSION` = 1,
 `ERR_MESSAGE` = 2,
 `ERR_MESSAGE_SIZE` = 3,
 `ERR_POINTER` = 4,
 `ERR_ACCESS` = 5,
 `ERR_ARGUMENT` = 6,
 `ERR_SIZE` = 7,
 `ERR_NOT_SUPPORTED` = 8,
 `ERR_TIMEOUT` = 9,
 `ERR_CHECKSUM` = 10,
 `ERR_OVERFLOW` = 11,
 `ERR_EVENT` = 12,
 `ERR_INTERRUPT` = 13,
 `ERR_TIMER_ACCESS` = 14,
 `ERR_LED_ACCESS` = 15,
 `ERR_TEMP_SENSOR_ACCESS` = 16,
 `ERR_DATA_TRANSFER` = 17,
 `ERR_FIFO` = 18,
 `ERR_OVER_TEMP` = 19,
 `ERR_IDENTIFICATION` = 20,
 `ERR_COM_INTERFACE` = 21,
 `ERR_SYNCHRONISATION` = 22,
 `ERR_PROTOCOL` = 23,
 `ERR_MEMORY` = 24,
 `ERR_THREAD` = 25,
 `ERR_SPI` = 26,
 `ERR_DAC_ACCESS` = 27,
 `ERR_I2C` = 28,
 `ERR_NO_DATA` = 29,
 `ERR_SYSTEM_CONFIG` = 30,
 `ERR_USB_ACCESS` = 31,
 `ERR_ADC_ACCESS` = 32,
 `ERR_SENSOR_CONFIG` = 33,
 `ERR SATURATION` = 34,
 `ERR_MUTEX` = 35,
 `ERR_ACCELEROMETER` = 36,
 `ERR_CONFIG` = 37,
 `ERR_BLE` = 38,
 `ERR_FILE` = 39,
 `ERR_DATA` = 40,
 `ERR_BUSY` = 41 }

2.3.1 Detailed Description

Generic error codes used by ams libraries.

2.3.2 Typedef Documentation

2.3.2.1 `err_code_t` typedef enum `error_codes` `err_code_t`

This definition will be used for function return values.

2.3.3 Enumeration Type Documentation

2.3.3.1 `error_codes` enum `error_codes`

Error Codes.

Enumerator

<code>ERR_SUCCESS</code>	Operation was successful.
<code>ERR_PERMISSION</code>	Operation is not permitted.
<code>ERR_MESSAGE</code>	Message is invalid. (Unsupported message type, incorrect CRC, ...)
<code>ERR_MESSAGE_SIZE</code>	Message has the wrong size.
<code>ERR_POINTER</code>	Pointer is invalid. (NULL pointer, pointer wrong memory region, ...)
<code>ERR_ACCESS</code>	Access is denied.
<code>ERR_ARGUMENT</code>	Argument is invalid.
<code>ERR_SIZE</code>	An argument has the wrong size.
<code>ERR_NOT_SUPPORTED</code>	Function is not supported or not implemented.
<code>ERR_TIMEOUT</code>	Operation timed out.
<code>ERR_CHECKSUM</code>	Checksum comparison failed.
<code>ERR_OVERFLOW</code>	Data overflow occurred.
<code>ERR_EVENT</code>	Getting or setting an event failed. (Event queue is full or empty, unexpected event received, ...)
<code>ERR_INTERRUPT</code>	Getting or setting an interrupt failed. (Interrupt resource is not available, ...)
<code>ERR_TIMER_ACCESS</code>	Accessing the timer peripheral failed.
<code>ERR_LED_ACCESS</code>	Accessing the LED peripheral failed.
<code>ERR_TEMP_SENSOR_ACCESS</code>	Accessing the temperature sensor failed.
<code>ERR_DATA_TRANSFER</code>	Communication error occurred.
<code>ERR_FIFO</code>	FIFO operation failed.
<code>ERR_OVER_TEMP</code>	Overtemperature detected.

Enumerator

ERR_IDENTIFICATION	Sensor identification failed.
ERR_COM_INTERFACE	Generic communication interface error. (Communication interface is not available, error while opening or closing a communication interface, ...)
ERR_SYNCHRONISATION	Synchronization error occurred.
ERR_PROTOCOL	Generic protocol error occurred.
ERR_MEMORY	Memory allocation error occurred.
ERR_THREAD	Thread handling operation failed.
ERR_SPI	Accessing the SPI peripheral failed.
ERR_DAC_ACCESS	Accessing the DAC peripheral failed.
ERR_I2C	Accessing the I2C peripheral failed.
ERR_NO_DATA	No data available.
ERR_SYSTEM_CONFIG	System configuration failed. (System resource is not available, system resource generates an error, ...)
ERR_USB_ACCESS	Accessing the USB peripheral failed.
ERR_ADC_ACCESS	Accessing the ADC peripheral failed.
ERR_SENSOR_CONFIG	Sensor configuration failed.
ERR_SATURATION	Saturation detected.
ERR_MUTEX	Mutex handling operation failed.
ERR_ACCELEROMETER	Accessing the accelerometer failed.
ERR_CONFIG	Software component is unusable due to incomplete or incorrect configuration.
ERR_BLE	BLE stack handling operation failed.
ERR_FILE	File handling operation failed.
ERR_DATA	Internal data inconsistency detected.
ERR_BUSY	Module is busy.

3 Data Structure Documentation

3.1 agc_configuration_t Struct Reference

Data Fields

- [agc_mode_t](#) mode
- [agc_ampl_cntl_mode_t](#) led_control_mode
- [agc_channel_id_t](#) channel
- [uint8_t](#) led_current_min
- [uint8_t](#) led_current_max
- [uint8_t](#) rel_amplitude_min_x100
- [uint8_t](#) rel_amplitude_max_x100
- [uint8_t](#) rel_amplitude_motion_x100
- [uint8_t](#) num_led_steps
- [uint8_t](#) reserved [3]
- [int32_t](#) threshold_min
- [int32_t](#) threshold_max

3.1.1 Detailed Description

Configuration of an instance of the AGC algorithm.

Note

For mode [AGC_MODE_STATUS](#) all the parameters but channel are unused and should be set to zero.

3.1.2 Field Documentation

3.1.2.1 mode [agc_mode_t](#) agc_configuration_t::mode

Selects the AGC algorithm mode.

3.1.2.2 led_control_mode [agc_ampl_cntl_mode_t](#) agc_configuration_t::led_control_mode

Selects the LED amplitude control mode.

3.1.2.3 channel [agc_channel_id_t](#) agc_configuration_t::channel

Selects the PPG channel that is controlled.

3.1.2.4 **led_current_min** `uint8_t agc_configuration_t::led_current_min`

Lower bound of the LED current range as a register value.

3.1.2.5 **led_current_max** `uint8_t agc_configuration_t::led_current_max`

Upper bound of the LED current range as a register value.

3.1.2.6 **rel_amplitude_min_x100** `uint8_t agc_configuration_t::rel_amplitude_min_x100`

Lower bound of the targeted PPG signal amplitude, relative to the size of the target PPG signal range. The unit of this field is percent. The minimum valid value is 0, the maximum valid value is 100.

3.1.2.7 **rel_amplitude_max_x100** `uint8_t agc_configuration_t::rel_amplitude_max_x100`

Upper bound of the targeted PPG signal amplitude, relative to the size of the target PPG signal range. The unit of this field is percent. This value must not be less than `rel_amplitude_min_x100`. The maximum valid value is 100.

3.1.2.8 **rel_amplitude_motion_x100** `uint8_t agc_configuration_t::rel_amplitude_motion_x100`

Minimum PPG signal amplitude at which the PPG signal is considered a motion artifact. The threshold is relative to the size of the target PPG signal range. The unit of this field is percent. This value must not be less than `rel_amplitude_max_x100`. The maximum valid value is 100.

3.1.2.9 **num_led_steps** `uint8_t agc_configuration_t::num_led_steps`

Number of steps the LED current range is partitioned in. When the AGC algorithm determines that the LED current needs to be increased or decreased and the bounds of the LED current range are not yet reached, the LED current is adjusted by one step. If `led_control_mode` is set to [AGC_AMPL_CNTL_MODE_AUTO](#), this value must not be zero and must not be greater than the size of the LED current range.

3.1.2.10 **reserved** `uint8_t agc_configuration_t::reserved[3]`

Padding bytes. Needs to be set to zero.

3.1.2.11 **threshold_min** `int32_t agc_configuration_t::threshold_min`

Lower bound of the target PPG signal range. The unit of this field is ADC counts.

3.1.2.12 **threshold_max** `int32_t agc_configuration_t::threshold_max`

Upper bound of the target PPG signal range. The unit of this field is ADC counts. This value must be greater than `threshold_min`.

3.2 agc_status_t Struct Reference

Data Fields

- [agc_change_state_t](#) `pd_offset_change`
- `uint8_t` `pd_offset_current`
- [agc_change_state_t](#) `led_change`
- `uint8_t` `led_current`

3.2.1 Detailed Description

Status information of an instance of the AGC algorithm.

3.2.2 Field Documentation

3.2.2.1 `pd_offset_change` [agc_change_state_t](#) `agc_status_t::pd_offset_change`

Change information for the photodiode offset current.

3.2.2.2 `pd_offset_current` `uint8_t` `agc_status_t::pd_offset_current`

Current value of the photodiode offset current as a register value.

3.2.2.3 `led_change` [agc_change_state_t](#) `agc_status_t::led_change`

Change information for the LED current.

3.2.2.4 `led_current` `uint8_t` `agc_status_t::led_current`

Current value of the LED current as a register value.

3.3 agc_version_t Struct Reference

Data Fields

- `uint8_t` `major`
- `uint8_t` `minor`
- `uint8_t` `patch`

3.3.1 Detailed Description

Version information of the AGC module.

3.3.2 Field Documentation

3.3.2.1 major `uint8_t agc_version_t::major`

Major version.

3.3.2.2 minor `uint8_t agc_version_t::minor`

Minor version.

3.3.2.3 patch `uint8_t agc_version_t::patch`

Patch version.